

SCENAR COSMODIC EX735Ag Modific

Professional Devices



SCENAR COSMODIC EX735Ag modific — is a professional SCENAR-COSMODIC series model, employing two technologies of treatment — SCENAR and COSMODIC in one device. Intended for both professional medical practice and home use. Fully automated action and operation.

“LET Medical” presents EX735Ag which sets the new standards of efficiency for professional models: two independent technologies in one device united by automatic selection system. Two healing technologies, SCENAR and COSMODIC, ensure that your approach to treatment is appropriate for each individual case and significantly reduce time needed for recovery. Presented with two models: EX735Ag AND EX735Ag Modific

Innovations

- automatically select healing technology which helps you to achieve optimal effect of treatment
- form subsequent impulses in advance, speeding up the recovery
- automatically regulate action energy levels in both SCENAR and COSMODIC which excludes “defense reaction effect” of the body
- correct manual settings of action modes if they are wrong
- start treatment immediately employing all functions of automatic control (switches on itself, tunes itself, etc)
- maintain perfect fine-tuning via user mode for automatic correction
- keep accurate records and regulate time of use of the device by medical personnel

SCENAR COSMODIC EX735Ag modific inherited the platform of its predecessor VX735Ag — the time-proved favourite allowing you to set SCENAR and COSMODIC healing technologies manually within one output stage. Now that these two technologies are separated into two independent output stages, the prioritized action is set automatically, providing deeper and more accurate treatment.

Electrodes

The electrodes are detachable and made in a form of three bars, which are flat on the underside and rounded at the front end of the device. The rounded part is convenient for treatment of small or hard-to-reach areas (for example, on the face). The flat part is mainly used for treatment of large areas. The electrodes are made of silver according to the special technology, which allows maintaining the natural structure of silver. This provides additional healing effect and more accurate signal of action.

Silver is known for its bactericidal, antiseptic, anti-inflammatory and styptic effects. Silver is a natural bactericidal metal effective against 650 kinds of bacteria which do not acquire resistance to it, like to most antibiotics. In addition, you can order electrodes made of special alloy of gold. Gold electrodes have the same shape as silver ones and can be easily installed in the device by the user.

Main features:

- touch switch-on;
- setting the preferred technology of treatment (SCENAR or COSMODIC) automatically;
- setting technology of treatment (SCENAR or COSMODIC) manually;
- determining areas for action (graphic and numerical gradation);
- regulating action energy automatically and manually;
- two action modes: treatment and cosmetology;
- “dynamic adaptation” mode allowing you to manage stimulation during prolonged action automatically;
- timer of action upon the area, timer of the whole procedure, timer of the resource of the device;
- switch-off timer;
- graphic OLED display;
- menu language: Russian, English
- LED and sound indication;
- switching sound off and on;
- vibrosignal;
- detachable silver electrodes;
- connecting remote electrodes;
- automatic tuning in to remote electrodes;
- automatic switch-over to action on mucous membranes when using rectal and vaginal electrodes;



Casing

The device is encased in “SMART” series casing made of shockproof (3.5 mm thick) ABS plastic. Waterproof keyboard.

Indication

OLED display

Graphic OLED display shows main action modes, action energy, activity of the area being acted upon, dosage, time of action upon the certain area and full time of action during session, menu and battery status.

LED indication

Additional LED indication allows you to evaluate if the area you have chosen for action is the appropriate one, activation of dynamic adaptation block and proper contact of the electrodes with skin.

Size

143 x 43 x 30 mm

Power supply

Two galvanic cells or rechargeable batteries, size AA, 1.5V.

Delivery set

- Device
- User’s Guide
- Package

Touch of Power Therapies

PO Box 239

Mapleton

Queensland 4560

Australia

www.touchofpower.com.au

ph 07 5445 7778

m 0407 136 794

For the Professional, this device is very easy and simple to operate. It is the best choice for home user if the price is not an issue.