

# SCENAR COSMODIC PS705

## Semi-professional Devices



**SCENAR COSMODIC PS705** is a semi-professional SCENAR COSMODIC series model employing COSMODIC technology. It is intended for both use at home and professional practice. Fully automated action and operation. Very power-saving.

### Innovations

The new model has all capabilities of pro devices. Full range of automated features — switch-on, regulation of energy and depth of action — always ensure that your treatment will be effective. Accurate indication of search of the optimal areas for treatment.

When you work with vaginal and rectal electrodes the device automatically switches to the special tender impulse and switches itself off after the minimum required time of action is over. This makes such delicate process simple and convenient in home environment.

Flat and thin casing of the device is convenient to lie on when you need prolonged action on the same area.

**SCENAR COSMODIC PS705 can be a wonderful assistant to the professional device.**

### Main features

- touch switch-on
- regulating action energy automatically
- “dynamic adaptation” mode allowing you to manage stimulation during prolonged action automatically
- determining areas for action
- LED and sound indication;
- switching sound off and on
- connecting remote electrodes;
- automatic tuning in to remote electrodes;

## Indication

### *LED indication displays:*

- contact of the electrodes with skin
- preliminary analysis of the treated area
- levels of action energy
- body reaction
- end of the action dose
- activation of “dynamic adaptation”
- battery status



## Electrodes

Electrodes of model PS705 are made of special steel which has a pyramidal structure. An important peculiarity of such kind of steel is that chrome does not penetrate into skin due to electrolysis process occurring during action.

Electrodes of model PS705 are placed on the front-end and lower part of the casing. The shape of the electrode allows placing the device in various planes which makes treatment of hard-to-reach areas easy.

The front-end part of the electrode is very convenient for treatment of small areas on the face and other awkward places. The flat part of the electrode on the lower side of the casing is convenient for treatment of large areas. The patient can lie on the electrode during prolonged action.

## Casing

SCENAR COSMODIC PS705Ag and PS705 are encased in “Pocket M” series casing. Waterproof keyboard and sealants make disinfection and cleaning of the device easy. The casing is made of shockproof (3 mm thick) ABS plastic.

## Size

115 x 58 x 23 mm

## Power supply

Two galvanic cells or rechargeable batteries, size AAA, 1.5V.

## Delivery set

- Device
- User’s Guide
- Package

The 705 series is the best choice for those who wish a robust device which would be simple to operate (tourism, extreme sports, first-aid kit in a car, patients to keep up between professional sessions, etc.). High quality and simple use are main highlights of this model

*Touch of Power Therapies*  
*PO Box 239*  
*Mapleton*  
*Queensland 4560*  
*Australia*  
*[www.touchofpower.com.au](http://www.touchofpower.com.au)*  
*ph 07 5445 7778*  
*m 0407 136 794*