

SCENAR COSMODIC PS705Ag mini

SCENAR COSMODIC PS705Ag mini is a personal-use SCENAR-COSMODIC series model employing COSMODIC technology. Fully automated action and operation. Very power-saving. The “smallest” COSMODIC. This is the case when size doesn’t matter — it is as capable as the “bigger” models.

Innovations

SCENAR COSMODIC PS705Ag mini developed on a platform of model PS705 mini. The main difference is silver electrodes with dotted magnets that increase the depth of influence. Algorithms of influence are also adjusted to functioning with silver electrodes, and are able to create a virtual wave pyramid with its vertex being inside of a body.

Main features

- touch switch-on
- regulating action energy automatically
- “dynamic adaptation” mode allowing you to manage stimulation during prolonged action automatically
- determining areas for action
- LED and sound indication;
- switching sound off and on

Indication

LED indication displays:

- contact of the electrodes with skin
- preliminary analysis of the treated area
- levels of action energy
- body reaction
- end of the action dose
- battery status



Electrodes

The silver electrodes of SCENAR COSMODIC PS705Ag mini are detachable and made in the form of three "buttons" on the back side of the case with dotted magnets that increase the depth of influence. The electrodes are made of silver according to the special technology, which allows maintaining the natural structure of silver. This provides additional healing effect and more accurate signal of action. Silver is known for its bactericidal, antiseptic, anti-inflammatory and styptic effects. Silver is a natural bactericidal metal effective against 650 kinds of bacteria which do not acquire resistance to it, like to most antibiotics.

Casing

SCENAR COSMODIC PS705Ag mini is encased in «Pocket S» series casing. Waterproof keyboard and sealants make disinfection and cleaning of the device easy. The casing is made of shockproof (3 mm thick) ABS plastic.

Size

85 x 46 x 16 mm

Power supply

Two galvanic cells or rechargeable batteries, size N, 1.5V.

Delivery set

- Device
- User’s Guide
- Package

Ideal for travelers and those who need good medical care away from home, the PS705mini is ready for every-day emergencies.

Touch of Power Therapies

PO Box 239

Mapleton

Queensland 4560

Australia

www.touchofpower.com.au

ph 07 5445 7778

m 0407 136 794